## Medical Solutions Customer Services

Name Gebhardt  
Dept. CS PS 1, (Med/LINA)  
Phone +49 9131 84 3319  
Fax +49 9131 84 7363  
E-Mail: MED CS PRODUCT TRACKING  
CS.Product\_Tracking@med.siemens.de  
Date: November 19, 2002

## Safety Modification SP018/02/S MODULARIS Uro Advisory Letter

### Reason for Update:

This customer safety advisory is to inform you that the **maximum weight capacity** of the **radiotransparent examination area** (refer to figure below) for the MODULARIS Uro is limited to **40 kg (90 lbs)** until new replacement parts are installed on your unit. Weights greater than the recommended 40 kg (90 lbs) on the examination area may cause damage to this part of the table, resulting in personal injury and, in some cases, serious bodily harm.

### Measures:

Immediately distribute the enclosed “**Letter to MODULARIS Uro Customers**” by registered mail to concerned customers and report distribution via the Siemens Intranet reporting tool by your authorized person.

http: [http://tdi-s3.med.siemens.de/ApplicationDB/SilverStream/Pages/up\\_tdi\\_ps\\_portal.html](http://tdi-s3.med.siemens.de/ApplicationDB/SilverStream/Pages/up_tdi_ps_portal.html)

### Systems affected:

Product:	Part No.	Serial No.
<b>MODULARIS Uro</b>	<b>55 31 012</b>	<b>01127 – 01155</b>

### Material:

No material is required. The referenced “**Letter to MODULARIS Uro Customers**” is enclosed with this letter.

### Deadline:

The deadline is February 19, 2003.

Very truly yours,

Mailing Address:  
Siemens AG  
CSPS 1 (Med/LINA)  
P.O.Box 3260  
91052 Erlangen

Geist

Gebhardt

Enclosure

To:

**Regional managers and CEO's of subsidiaries**

Name	Markus Lanski
Department	SPU
Telephone	+49 9131 84-6644
Fax	+49 9131 84-5400
E-Mail	markus.lanski@siemens.com
Our reference	
Date	November 28, 2002

### **URGENT! For immediate distribution!**

#### **Customer safety advisory re: patient positioning on Siemens MODULARIS Uro systems**

Dear Sir or Madam,


Please forward the enclosed advisory letter immediately to the MODULARIS Uro customers in your region to whom this pertains. After reading the letter, file it in the safety register of the operating manual for the respective customer(s). Refer to the attachment for a list of the MODULARIS Uro systems affected.

As a result of faulty table insert parts (i.e. the layer in the 'sandwich' construction of the radiotransparent examination area) we are obligated to replace these parts on the respective systems in operation in your sales region.

For additional details, refer to the enclosed sample letter. (For legal reasons, the letter may not be modified.) You are responsible for translating this customer letter into the national language and informing customers without delay. Please confirm that you have received this letter by returning a copy to our offices.

We wish to thank you in advance for your understanding and support in updating our MODULARIS Uro customers.

With best regards

Markus Lanski  
Product Manager for Lithotripsy

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To:

**All MODULARIS Uro customers affected**

Name	Markus Lanski
Department	SPU
Telephone	+49 9131 84-6644
Fax	+49 9131 84-5400
E-Mail	markus.lanski@siemens.com
Our reference	
Date	November 28, 2002

## Customer Safety Advisory


**Re: MODULARIS Uro patient table**

**Weight capacity of the radiotransparent area limited to 40 kg (90 lbs) max.**

Dear MODULARIS Uro customer,

This customer safety advisory is to inform you that the **maximum weight capacity** of the **radiotransparent examination area** (refer to figure below) for the MODULARIS Uro is limited to **40 kg (90 lbs)** until new replacement parts are installed on your unit.

Weights greater than the recommended 40 kg (90 lbs) on the examination area may cause damage to this part of the table, resulting in personal injury and, in some cases, serious bodily harm.

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In addition, please review the safety notice in the operating manual (Safety register, section 2.11) concerning the proper access position (refer to figure). Patient movements should be limited on the transfer area due to hazards mentioned in the manual.

The maximum patient weight for the entire table is 180 kg (400 lbs).

Siemens Service will replace the insert parts as soon as possible and at no cost to you.  
Until then, please use increased caution when operating your MODULARIS Uro.

Please inform your staff about the contents of this customer safety advisory and instruct them to comply with its guidelines.

Thank you for your understanding and cooperation in this matter.

With best regards

---

Markus Lanski  
Product Manager for Lithotripsy

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